

PeliPIP Instructions for use

Choose the composition sheet you wish.

Look for cell suspension and batch number:



0,5-11P_8000080621-SCR_8000080619-16P_8000080621.xml

0,5: 0.5% cell suspension

11P: 11 cells panel; these are the same as the first 11 cells of the 16 cells panel

80000...: batch number

SCR: screening panel

16P: 16 cells panel

- Click with the right mouse button on the desired composition sheet / file.
- Click on "save target as"
- Save the chosen composition sheet / file in a folder of your choice.

Do not change the file name!

- Click on "close".

You can upload this file in PeliPIP